
**MIKRALGEN**

"La Nitrière" - 01160 Neuville sur Ain

FRANCE

Tel.: 33 (0)4 74 37 72 28

contact@mikralgen.com

English

# KIT HYGI-PRO

Italian

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*Colorimetric determination of proteins and reducing oses*

*Qualitative control of surface hygiene and rinse water*

Hygi-pro kit is used for detecting residuals of proteins and reducing oses on surfaces. This colorimetric determination indicates the hygiene level of surfaces, this technique is easy to carry out and doesn't require a sophisticated equipment. Packaging tubes improve efficiency and save material.

## PRINCIPLE

In alkaline conditions, proteins and reducing oses reduce cupric ions ( $\text{Cu}^{++}$ ) in ( $\text{Cu}^+$ ) forming complex with the modified Biuret reagent later revealed by coloured reaction.

Thus it is possible to determine the protein quantity and reducing oses by examining tone of the color.

**KIT : 28 swabs, Reagent A ready to use, 24 coated tubes**

## STABILITY / SENSITIVITY

Store at room temperature, not in contact with sun-light, in a closed box. Detection limit of 50-100µg proteins and/or reducing oses per swabing.

- **Do not touch materials with fingers because many proteins contaminate skin** and that would interfere with reaction developed in Hygi-pro test.
- **Keep away reagents from sun-light.**
- **Use reagents with care to not contaminate them.**
- **Irritant product for skin and eyes** ; in case of skin or eyes or mucous membrane contact with liquids, rinse immediately with a lot of water.

## OPERATING INSTRUCTIONS

- 1/ ♦ Open the aluminium package to take number tubes you need, you can separate the strip of plugs with scissors.
  - ♦ Identify the test sample by writing on the top of the tube, fit the required number of coated tubes onto the holder ; we recommend to use one tube as 'negative control' for each test sample series.
  - ♦ Fill the required number of coated tube with reagent A until the small external relief on the tube. **Keep one tube as the negative control.** The prepared solution in the reconstituted tubes must be used within 20 minutes maximum.
- 2/ **Sample preparation : surface - rinse water**
  - ♦ **Surface** : If the surface is dry, we recommend to humidify the surface or directly the sampling head of the swab with one drop of reagent A.
  - Swab the surface you want to test by turning up the swab, an area of 10 cm x 10 cm is recommended in food processing industries. Carry out the manipulation by reading step 3.
  - ♦ **Rinse water sample** : Put the swab in the water you want to test.
- 3/ **Test procedure**
  - ♦ Put the sample swab into the sample tube prepared in step 1.
  - ♦ **Homogenize gently the solution by stirring the swab in the tube. Incubate at room temperature (20-25°C) for 10 minutes. During the incubation step, let the swab in the sample tube on the holder.**
- 4/ **Reaction procedure**
  - Incubate for 10 min.** at room temperature (20-25°C).

## INTERPRETATION / RESULTS

- ♦ Read the result after a 10 minutes incubation period. **Homogenize gently solution by turning the swab in the tube before interpretation.** The colour of the negative control tube must be light-green and the positive test has turned purple with variable tone. The more the purple is dark, the more the sample contains proteins and/or reducing oses. Control the light-green colour of the negative tube to validate the manipulation.
  - ♦ Compare the tests tubes with the color-checker supplied to evaluate residuals proteins and/or reducing oses detected on the surfaces tested.
  - ♦ The colouration intensity of positive tests increases with time reaction and temperature. It could be possible to detect protein concentration less than 50µg per sample by increasing the incubation period to 20 or 30 minutes at 37°C.
- All Hygi-pro tests must be conducted at the same working conditions (incubation period, temperature) to carry out reproducibility in order to allowing conclusions in hygiene monitoring procedures.**